

# Anti-LINE-1 ORF1p Rabbit antibody



Catalog Number: 167990

## Product name

Anti-LINE-1 ORF1p Rabbit antibody

## Specificity

Human

## Antibody description

LINE-1 ORF1p Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human LINE-1 ORF1p

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/50

IP: 1/20

## Validations



Western blot detection of LINE-1 ORF1p in HeLa,A549,HL-60,U2OS cell lysates using LINE-1 ORF1p Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:42kDa.