

Anti-eIF4EBP1 Rabbit antibody



Catalog Number: 168990

Product name

Anti-eIF4EBP1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

eIF4EBP1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human eIF4EBP1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/5000


IHC: 1/50

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of eIF4EBP1 in Jurkat, C6, CHO-K1 cell lysates using eIF4EBP1 Rabbit mAb(1:1000 diluted). Predicted band size:13kDa.Observed band size:15-20kDa.