

ECD antibody

Catalog Number: 110290

Product name

ECD antibody

IP: 1:200-1:1000

IF: 1:10-1:100


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

ECD Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 73 kDa




Preparation

This antibody was obtained by immunization of ECD recombinant protein (Accession Number: NM_007265). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110290(ECD antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

IP Result of anti-ECD (IP:Catalog No:110290, 4ug; Detection:Catalog No:110290 1:300) with HeLa cells lysate 2480ug.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunofluorescent analysis of HeLa cells using Catalog No:110290(ECD Antibody) at dilution of

ECD antibody



Catalog Number: 110290

1:25 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)