

# ARL6IP5 Polyclonal Antibody



Catalog Number: 162290

## Product name

ARL6IP5 Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to ARL6IP5

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human ARL6IP5 (NP\_006398.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - ARL6IP5 Polyclonal Antibody

Western blot analysis of extracts of mouse thymus, using ARL6IP5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.