

# LNX1 antibody



Catalog Number: 112290

---

## Product name

LNX1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LNX1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse kidney tissue. Observed molecular weight by Western blot: 70 kDa

## Preparation

This antibody was obtained by immunization of LNX1 recombinant protein (Accession Number: XM\_017008776). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112290(LNX1 antibody) at dilution of 1:70