

# Rabbit anti DDDDK-Tag pAb



Catalog Number: 166290

## Product name

Rabbit anti DDDDK-Tag pAb

## Antibody description

Polyclonal antibody to Rabbit anti DDDDK-Tag

## Preparation

Antigen: Recombinant Protein of DDDDK

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF, IP

## Dilutions


WB 1:1000 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100


IP 1:50 - 1:100

## Validations



Western blot - Rabbit anti DDDDK-Tag pAb

Western blot analysis of recombinant DDDDK protein with Rabbit anti DDDDK-Tag pAb at dilution of 1:5000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - Rabbit anti DDDDK-Tag pAb

Immunofluorescence analysis of HeLa cells transfected with DDDDK-Tag (top left) and untreated HeLa cells (top right) use Rabbit anti DDDDK-Tag (red,) at dilution of 1:100. Blue: DAPI for nuclear staining.