

CHST1 antibody



Catalog Number: 109290

Product name

Polyclonal

CHST1 antibody

Ig Type

Specificity

Rabbit IgG

Human,Mouse,Rat; other species not tested.

Applications

Antibody description

ELISA, WB

CHST1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 48-55 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of CHST1 recombinant protein (Accession Number: NM_003654). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109290(CHST1 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality