

# IgG control antibody

Catalog Number: 111690

## Product name

IgG control antibody

## Antibody description

IgG control Rabbit Polyclonal antibody. Positive IP detected in Control. Positive WB detected in Protein, Control.

## Preparation

This antibody was obtained by immunization of N/A. Purification method: Protein A purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IP: N/A

## Validations



IP sample detected with different secondary antibodies. Normal Rabbit IgG (Catalog No:111690) as control.



IP sample detected with different secondary antibodies. Normal Rabbit IgG (Catalog No:111690) as control.