

# Anti-Phospho-Nucleolin (Thr84) Rabbit antibody



Catalog Number: 168690

## Product name

Anti-Phospho-Nucleolin (Thr84) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB

## Antibody description

Phospho-Nucleolin (Thr84) Rabbit monoclonal antibody


## Dilutions

WB: 1/1000-1/5000

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr84 of human Nucleolin

## Validations



## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

Western blot detection of Phospho-Nucleolin (Thr84) in U251MG, HepG2, Hela cell lysates using Phospho-Nucleolin (Thr84) Rabbit mAb(1:500 diluted). Predicted band size: 77kDa. Observed band size: 100kDa.