

PRKAG3 Polyclonal Antibody



Catalog Number: 165390

Product name

PRKAG3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to PRKAG3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PRKAG3 (NP_059127.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:100

Validations



Immunohistochemistry - PRKAG3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using PRKAG3 antibody at dilution of 1:100 (40x lens).