

Anti-ZSCAN25 antibody



Catalog Number: 176390

Product name

Anti-ZSCAN25 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZSCAN25

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf498. aa range:251-300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of A549 cells using ZNF498 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).