

TRIP4 antibody



Catalog Number: 116390

Product name

TRIP4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRIP4 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, COLO 320 cells, HeLa cells. Observed molecular weight by Western-blot: 60-65kd

Preparation

This antibody was obtained by immunization of TRIP4 recombinant protein (Accession Number: NM_016213). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116390(TRIP4 antibody) at dilution of 1:1000