

# FUCA1 antibody

Catalog Number: 110790

## Product name

FUCA1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FUCA1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human heart tissue, human kidney tissue, human liver cancer tissue, human placenta tissue, human testis tissue.

## Preparation

This antibody was obtained by immunization of FUCA1 recombinant protein (Accession Number: BC017338). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:110790(FUCA1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110790(FUCA1 antibody) at dilution of 1:50 (under 40x lens)