

# ING1 antibody



Catalog Number: 111790

## Product name

ING1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ING1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A375 cells. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of ING1 recombinant protein (Accession Number: NM\_198219). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111790(ING1 Antibody) at dilution of 1:600