

# PGM2 antibody

Catalog Number: 113790

## Product name

PGM2 antibody

## Specificity

human; other species not tested.

## Antibody description

PGM2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HEK-293 cells. Positive IHC detected in human pancreas cancer. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of PGM2 recombinant protein (Accession Number: NM\_018290). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113790(PGM2 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113790(PGM2 antibody) at dilution of 1:100 (under 10x lens)