

Syntaxin 10 antibody



Catalog Number: 115790

Product name

Syntaxin 10 antibody

Immunogen

[Human STX10 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

Syntaxin 10 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of Syntaxin 10 recombinant protein (Accession Number: NM_001271611). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115790(STX10 antibody) at dilution of 1:200