

GNPDA1 antibody



Catalog Number: 111008

Product name

GNPDA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GNPDA1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in K-562 cells, HeLa cells, human brain tissue, rat kidney tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of GNPDA1 recombinant protein (Accession Number: XM_047416582). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111008(GNPDA1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111008(GNPDA1 antibody) at dilution of 1:50 (under 10x lens)