

# PLRG1 antibody

Catalog Number: 114008

## Product name

PLRG1 antibody

## Specificity

Human; other species not tested.

## Antibody description

PLRG1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in K-562 cells, BxPC-3 cells, HEK-293 cells. Observed molecular weight by Western-blot: 55-60kd

## Preparation

This antibody was obtained by immunization of PLRG1 recombinant protein (Accession Number: NM\_001201564). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114008(PLRG1 antibody) at dilution of 1:1000



IP Result of anti-PLRG1 (IP:Catalog No:114008, 4ug; Detection:Catalog No:114008 1:1000) with HEK-293 cells lysate 1400ug.