

ZNF593 antibody



Catalog Number: 117008

Product name

ZNF593 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF593 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, PC-3 cells. Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of ZNF593 recombinant protein (Accession Number: BC002580). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:117008(ZNF593 antibody) at dilution of 1:200