

Anti-BAT3/BAG-6 Rabbit antibody



Catalog Number: 169008

Product name

Rabbit IgG

Anti-BAT3/BAG-6 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000-1/5000

BAT3/BAG-6 Rabbit monoclonal antibody

ICC/IF: 1/50

Preparation

FC: 1/20

Antigen: A synthetic peptide of human BAT3

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of BAT3/BAG-6 in K562, Rat Brain, C6, 3T3 cell lysates using BAT3/BAG-6 Rabbit mAb(1:1000 diluted). Predicted band size: 119kDa. Observed band size: 150kDa.