

LYNX1 antibody

Catalog Number: 112408

Product name

LYNX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LYNX1 Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus cancer tissue, human oesophagus tissue, mouse brain tissue.

Preparation

This antibody was obtained by immunization of LYNX1 recombinant protein (Accession Number: BC032306). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human oesophagus cancer tissue slide using Catalog No:112408(LYNX1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human oesophagus cancer tissue slide using Catalog No:112408(LYNX1 Antibody) at dilution of 1:50 (under 40x lens)