

# Anti-DDX3 (6G8) Mouse antibody

Catalog Number: 168408

## Product name

Anti-DDX3 (6G8) Mouse antibody

## Specificity

Human, Rat, Mouse, Monkey

## Antibody description

DDX3 (6G8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human DDX3 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


WB, IP, ICC/IF, ChIP

## Dilutions

WB: 1/1000


;ICC: 1/200

## Validations




Western blot detection of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat cell lysate using DDX3 mouse mAb (1:1000 diluted). Predicted band size: 75KDa. Observed

band size: 75KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using DDX3



mouse mAb (dilution 1:200).

Immunoprecipitation analysis of HeLa cell lysates using DDX3 mouse mAb.