

RPIA antibody

Catalog Number: 114808

Product name

RPIA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPIA Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, HEK-293 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of RPIA recombinant protein (Accession Number: NM_144563). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:114808(RPIA antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using RPIA antibody Catalog No:114808 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).