

Anti-TRIM72 antibody

Catalog Number: 175808

Product name

Anti-TRIM72 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TRIM72

Preparation

This antigen of this antibody was recombinant protein of trim72

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:1000

IHC-P: 1:200-1:500

Validations




Fig1:: Western blot analysis of Mouse skeletal muscle using TRIM72 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded Mouse Pancreas Tissue using TRIM72 Polyclonal Antibody.