

Recombinant Human KRAS Protein



Catalog Number: 552108

Gene Name Synonym

PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A

Protein Construction

DNA encoding 1-169 aa (P01116-2) were fused with N-terminal 6His tag.

Organism

Human

Expression Host

E.coli

Activity

Not tested.

Purity

>90%, by SDS-PAGE with Coomassie Brilliant Blue staining.

Endotoxin

Please contact the lab for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

Predictes a molecular mass of 21.55 kDa. In SDS-PAGE under reducing conditions, it migrates as an

approximately 23.7 kDa band.

Formulation

Lyophilized from sterile PBS, pH 7.4


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

SDS-PAGE



Recombinant KRAS protein (His tag) SDS-PAGE

M: Marker; Lane 1: KRAS; Lane 2: BSA