

Anti-EGFR (1A11) Mouse antibody



Catalog Number: 168108

Product name

WB, IP

Anti-EGFR (1A11) Mouse antibody

Dilutions

Specificity

WB: 1/2000

Human

Antibody description

EGFR (1A11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human EGFR protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 0.1%BSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates



using EGFR mouse mAb.

Western blot detection of EGFR in A549 cell lysates using EGFR mouse mAb(dilution 1:2000).Predicted band size:134 Kda.Observed band size:175KDa.