

# Anti-Myc-tag (9C4) Mouse antibody



Catalog Number: 169108

## Product name

Anti-Myc-tag (9C4) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Tag

## Applications

WB, IP

## Antibody description

Myc-tag (9C4) Mouse monoclonal antibody


## Dilutions

WB:1/1000

## Preparation

Antigen: Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (EQKLISEEDL) coupled to KLH.

## Validations



Western blot analysis of extracts from CHO-K1 cells (A) or CHO-K1 cells transfected with different Myc-fusion proteins (Bu3001C), using Myc-tag mouse mAb (1:1000 diluted).

## Formulation

1×PBS with 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal