

MBNL3 antibody



Catalog Number: 112508

Product name

MBNL3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MBNL3 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of MBNL3 recombinant protein (Accession Number: NM_001170704). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112508(MBNL3 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:112508(MBNL3 Antibody) at dilution of 1:50 (under 40x lens)