

RNH1 antibody



Catalog Number: 107508

Product name

RNH1 antibody

Immunogen

[Human RNH1 Recombinant protein \(His tag\)](#)

Specificity

human; other species not tested.

Antibody description

RNH1 Mouse Monoclonal antibody. Positive WB detected in lo2 cells. Observed molecular weight by Western-blot: 50kD

Preparation

This antibody was obtained by immunization of RNH1 recombinant protein (Accession Number: NM_203388). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107508(RNH1 antibody) at dilution of 1:500