

CPI17 antibody



Catalog Number: 109508

Product name

CPI17 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CPI17 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive WB detected in HepG2 cells, HeLa cells, human colon tissue, PC-3 cells. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of CPI17 recombinant protein (Accession Number: NM_033256). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109508(PPP1R14A antibody) at dilution of



Immunofluorescent analysis of PC-3 cells using Catalog No:109508(PPP1R14A Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)