

# PIK3CA Polyclonal Antibody



Catalog Number: 164908

## Product name

PIK3CA Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PIK3CA

## Preparation

Antigen: Recombinant protein of human PIK3CA

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IP, FC

## Dilutions


WB 1:500 - 1:2000

IHC 1:20 - 1:200

IP 1:20 - 1:50

FC 1:20 - 1:50

## Validations



Western blot - PIK3CA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PIK3CA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.