

RPS4X antibody

Catalog Number: 114908

Product name

RPS4X antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

RPS4X Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells, human liver tissue, MCF7 cells, mouse pancreas tissue. Positive IHC detected in human pancreas cancer tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of RPS4X recombinant protein (Accession Number: NM_001007). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114908(SCAR antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114908(SCAR antibody) at dilution of 1:100 (under 10x lens)

RPS4X antibody


Catalog Number: 114908



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114908(SCAR antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells, using RPS4X antibody Catalog No:114908 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-SCAR (IP:Catalog No:114908, 3ug; Detection:Catalog No:114908 1:500) with HepG2 cells lysate 3200ug.