

ZBTB10 antibody



Catalog Number: 116908

Product name

ZBTB10 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZBTB10 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 65 kDa; 130 kDa

Preparation

This antibody was obtained by immunization of ZBTB10 recombinant protein (Accession Number: BC127098). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:500-1:5000

Validations



IP Result of anti-ZBTB10 (IP:Catalog No:116908, 3ug; Detection:Catalog No:116908 1:1000) with



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:116908(ZBTB10 Antibody) at dilution of 1:300