

FGF18 antibody



Catalog Number: 107208

Product name

FGF18 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FGF18 Mouse Monoclonal antibody. Positive WB detected in HT-1080 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of FGF18 recombinant protein (Accession Number: NM_033649). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:107208(FGF18 Antibody) at dilution of 1:1000