

# IDI2 antibody

Catalog Number: 111608

---

## Product name

Polyclonal

IDI2 antibody

## Ig Type

## Specificity

Rabbit IgG

human; other species not tested.

## Applications

## Antibody description

ELISA, WB

IDI2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 27-30kd

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Preparation

This antibody was obtained by immunization of IDI2 recombinant protein (Accession Number: NM\_033261). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111608(IDI2 antibody) at dilution of 1:2000

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality