

PAX7 antibody

Catalog Number: 113608

Product name

PAX7 antibody

Specificity

Human; other species not tested.

Antibody description

PAX7 Rabbit Polyclonal antibody. Positive WB detected in NCCIT cell, C2C12 cell, HeLa cells, Jurkat cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002584). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113608(PAX7 Antibody) at dilution of 1:1000



NCCIT cell were subjected to SDS PAGE followed by western blot with Catalog No:113608(PAX7 Antibody) at dilution of 1:600