

Anti-HNF-1-alpha/HNF1A antibody



Catalog Number: 104608

Product name

Anti-HNF-1-alpha/HNF1A antibody

Specificity

Human HNF-1-alpha / HNF1A

Antibody description

Rabbit polyclonal to HNF-1-alpha/HNF1A

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human HNF-1-alpha / HNF1A, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF, ICC/IF

Dilutions

ICC/IF: 0.5-1.5 µg/mL

Validations



HNF-1-alpha / HNF1A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence

Immunofluorescence staining of HNF1A in HeLa cells. Cells were fixed with 4% PFA, blocked with 10% serum. Then incubated with rabbit anti-human HNF1A polyclonal antibody (1 µg/ml)(Figure A), incubated with rabbit anti-human HNF1A polyclonal antibody and antigen (Figure B) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to nucleus.