

Anti-TrkA/NTRK1 antibody



Catalog Number: 106608

Product name

Anti-TrkA/NTRK1 antibody

Immunogen

[Rat TrkA/NTRK1 \(His Tag\) recombinant protein](#)

Specificity

Rat TrkA / NTRK1

Antibody description

Rabbit polyclonal to TrkA/NTRK1

Preparation

Produced in rabbits immunized with purified, recombinant Rat TrkA / NTRK1 (rh TrkA / NTRK1; NTRK1_RAT; Met1-Pro418). TrkA / NTRK1 specific IgG was purified by Rat TrkA / NTRK1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 µg/mL


This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat TrkA / NTRK1. The detection limit for Rat TrkA / NTRK1 is < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

Validations



TrkA / NTRK1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



TrkA / NTRK1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry