

USP34 antibody

Catalog Number: 116608

Product name

USP34 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

USP34 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human prostate cancer tissue. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, COLO 320 cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC022783). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116608(USP34 antibody) at dilution of 1:600



IP Result of anti-USP34 (IP:Catalog No:116608, 4ug; Detection:Catalog No:116608 1:500) with mouse brain tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

USP34 antibody

Catalog Number: 116608

No:116608(USP34 Antibody) at dilution of 1:50
(under 10x lens)



Immunohistochemistry of paraffin-embedded
human breast cancer tissue slide using Catalog
No:116608(USP34 Antibody) at dilution of 1:50
(under 40x lens)