

STIM1 antibody

Catalog Number: 107608

Product name

STIM1 antibody

Specificity

Rat, Human, Mouse; other species not tested.

Antibody description

STIM1 Mouse Monoclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells, K-562 cells, RAW 264.7 cells. Observed molecular weight by Western-blot: 77 kDa

Preparation

This antibody was obtained by immunization of STIM1 recombinant protein (Accession Number: NM_003156). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107608(STIM1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HEK-293 cells using Catalog No:107608(STIM1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)