

Anti-LAT Rabbit antibody



Catalog Number: 168608

Product name

Anti-LAT Rabbit antibody

Specificity

Human

Antibody description

LAT Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human LAT

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of LAT in U2OS cell lysates using LAT Rabbit mAb(1:500 diluted). Predicted band size:28kDa. Observed band size:36kDa.