

EEA1 antibody



Catalog Number: 110308

Product name

EEA1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EEA1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells, HepG2 cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003566).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Catalog No:110308(EEA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)