

HRP-conjugated Mouse anti His-Tag mAb



Catalog Number: 166308

Product name

HRP-conjugated Mouse anti His-Tag mAb

Applications

WB

Antibody description

Monoclonal antibody to HRP-conjugated Mouse anti His-Tag

Dilutions

WB 1:2000 - 1:10000


Preparation

Antigen: A synthetic peptide of human His-Tag

Validations

Formulation

PBS, 50% glycerol, pH7.3.



Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Western blot - HRP-conjugated Mouse anti His-Tag mAb

Ig Type

Mouse IgG

Western blot analysis of recombinant His protein with HRP-conjugated Mouse anti His-Tag mAb at dilution of 1:5000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.