

Anti-BMP4 Rabbit antibody



Catalog Number: 167308

Product name

Anti-BMP4 Rabbit antibody

WB: 1/500-1/2000

IHC: 1/50-1/200

Specificity

Human

ICC/IF: 1/50-1/200

IP: 1/50

FC: 1/50

Antibody description

BMP4 Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthesized peptide derived from human BMP4

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP, FC

Dilutions



Western blot analysis of extracts from THP-1 and Raji cells using BMP4 Rabbit pAb at 1:1000 dilution. Predicted band size: 47kDa. Observed band size: 47kDa.



Immunohistochemical analysis of paraffin-embedded human breast cancer, using BMP4 Antibody.