

HRP-conjugated GST Tag antibody



Catalog Number: 107308

Product name

Monoclonal

HRP-conjugated GST Tag antibody

Ig Type

Specificity

Mouse IgG2a

Tag

Applications

Antibody description

ELISA, WB

HRP-conjugated GST Tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: -

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: FN315687). Purification method: Protein A Purified.

Validations



Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Recombinant protein were subjected to SDS PAGE followed by western blot with HRP-66001(GST Tag Antibody) at dilution of 1:10000

Clonality