

HTR1A antibody



Catalog Number: 111708

Product name

HTR1A antibody

Specificity

Human; other species not tested.

Antibody description

HTR1A Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000524). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111708(HTR1A Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111708(HTR1A Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).