

# ACTL7A antibody



Catalog Number: 107708

## Product name

ACTL7A antibody

## Immunogen

[Human ACTL7A Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ACTL7A Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of ACTL7A recombinant protein (Accession Number: NM\_006687). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107708(ACTL7A antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using ACTL7A antibody Catalog No:107708 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).