

C21orf66 antibody



Catalog Number: 108708

Product name

C21orf66 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C21orf66 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 93 kDa

Preparation

This antibody was obtained by immunization of C21orf66 recombinant protein (Accession Number: BC030539). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108708(C21orf66 antibody) at dilution of



Immunofluorescent analysis of Hela cells, using C21orf66 antibody Catalog No:108708 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).