

Drebrin antibody



Catalog Number: 110018

Product name

Drebrin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Drebrin Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, mouse brain tissue, Neuro-2a cells, SH-SY5Y cells. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 100-130 kDa

Preparation

This antibody was obtained by immunization of Drebrin recombinant protein (Accession Number: NM_004395). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110018(Drebrin Antibody) at dilution of 1:3000



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:110018(Drebrin Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)