

# Anti-PGBD5 antibody



Catalog Number: 175018

## Product name

WB, ICC

Anti-PGBD5 antibody

## Dilutions

## Specificity

WB: 1:10,000–1:20,000

Human

ICC: 1:100–1:200

## Antibody description

Rabbit polyclonal antibody to PGBD5

## Preparation

This antigen of this antibody was synthetic peptide within human piggybac transposable element-derived protein 5 aa 61-115.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

## Validations




Fig1: Western blot analysis on D3 using anti-PGBD5 polyclonal antibody.




Fig2: ICC on HepG2 using anti-PGBD5 polyclonal antibody.